






## Product data sheet REF 57.36.25

11.12.2025

General data			
article number	57.36.25	Toennis dissector d/e 3.5/5mm 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11290	GTIN	04046826083458
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;04;aj;5X
Net weight (kg)	0,020	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
			
size	3.5/5mm 24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		