








Product data sheet REF 57.36.51

11.12.2025

General data			
article number	57.36.51	Penfield dissector d/e fig 1 17.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11290	GTIN	04046826083489
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;04;aj;5X
Net weight (kg)	0,019	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	17.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		