



## Product data sheet REF 57.37.02

11.12.2025

General data			
article number	57.37.02	MC Micro Cushing rong 17.0cm DieBlauen.2 2x5mm str serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826449414
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.04.58 (Medicon DieBLAUEN.2 dism Rongeurs w Silicon Handle)		
Complete picture		Detail images	
 			
size	17.0cm 2x5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	ELASTOSIL LR3003 (Silikonkautschuk)		