



Product data sheet REF 57.37.04

11.12.2025

| General data | | | |
|---|--|---|---|
| article number | 57.37.04 | MC Cushing rong DieBlauen.2 18cm 2x10mm 30° upw | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35541 | GTIN | 04046826449438 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-669 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;06;ai;5X |
| Net weight (kg) | 0,000 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-2 |
|   | | | |
| Medicon media | 451.04.58 (Medicon DieBLAUEN.2 dism Rongeurs w Silicon Handle) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 18cm 2x10mm 30° | Coating | Artificially created PVD layer of titanium aluminum nitrite e.g. Magic Cut scissors |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Handle | ELASTOSIL LR3003 (Silikonkautschuk) | | |