







Product data sheet REF 57.37.05

11.12.2025

General data			
article number	57.37.05	Micro Cushing rong DieBlauen.2 17.0cm 2x5mm 30° ang upw serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826467982
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.04.58 (Medicon DieBLAUEN.2 dism Rongeurs w Silicon Handle)		
Complete picture		Detail images	
 			
size	17.0cm 2x5mm 30°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	ELASTOSIL LR3003 (Silikonkautschuk)		