



Product data sheet REF 57.37.57

11.12.2025

General data			
article number	57.37.57	Caspar rong DieBlauen.2 15.5cm 4x14mm str fenest serr	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826449643
Risk category	1r (EU MDR)	UMDNS-Code	15-669
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-2
			
Medicon media	451.04.58 (Medicon DieBLAUEN.2 dism Rongeurs w Silicon Handle)		
Complete picture		Detail images	
			
size	15.5cm 4x14mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw Handle	1.4021 (X20Cr13)		
	1.4305 (X8CrNiS18-9)		
	ELASTOSIL LR3003 (Silikonkautschuk)		