








## Product data sheet REF 57.49.10

11.12.2025

General data			
article number	57.49.10	Eif hook 90° sh 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	46768	GTIN	04046826084677
Risk category	1r (EU MDR)	UMDNS-Code	12-029
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,100	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.53 (catalogue cardio )		
Complete picture		Detail images	
  			
size	90° 24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		