








Product data sheet REF 57.52.05

11.12.2025

General data			
article number	57.52.05	Cushing retr 20cm 10mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826085049
Risk category	1r (EU MDR)	UMDNS-Code	12-029
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,068	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	20cm 10mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		