





Product data sheet REF 57.55.37

11.12.2025

General data			
article number	57.55.37	Williams retractor 17cm 50x20mm rt/hook lf	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826421588
Risk category	2a (EU MDD)	UMDNS-Code	13-384
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media			
Complete picture		Detail images	
			
size	17cm 50x20mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		