





## Product data sheet REF 57.56.41

11.12.2025

| General data   |                                   |                                  |                |
|--|-----------------------------------|----------------------------------|----------------|
| article number   | 57.56.41                          | Cloward blade toothed 40x15/23mm |                |
| Unit   | piece                             | Customs code (HS code)           | 90189084       |
| GMDN   | 45182                             | GTIN                             | 04046826086008 |
| Risk category  | 2a (EU MDD)                       | UMDNS-Code                       | 13-390         |
| Classification   | 7 (EU MDD)                        | Base UDI-DI                      |                |
| Net weight (kg)  | 0,015                             | Operating instructions           |                |
| Latex free   | yes                               | Preparation instructions         | MPK01          |
|   |                                   |                                  |                |
| Medicon media  |                                   |                                  |                |
| Complete picture   |                                   | Detail images                    |                |
|  |                                   |                                  |                |
| size   | 40x15/23mm                        | Coating                          |                |
| Storage conditions   |                                   |                                  |                |
| No special requirements  |                                   |                                  |                |
| Materials  |                                   |                                  |                |
|  | Material number (DIN designation) |                                  |                |
| Main part  | 1.4301 (X5CrNi18-10)              |                                  |                |
| Others   | 1.4305 (X8CrNiS18-9)              |                                  |                |