






Product data sheet REF 57.58.35

11.12.2025

| General data | | | |
|--|-----------------------------------|--------------------------|--|
| article number | 57.58.35 | Piccolino speculum 40mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45182 | GTIN | 04046826383961 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 13-390 |
| Classification | 7 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,000 | Operating instructions | G575752 |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 40mm | Coating | Artificially created PVD layer of titanium nitrite e.g. Blackline scissors |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |