







Product data sheet REF 57.64.40

11.12.2025

General data			
article number	57.64.40	Gigli (Olivecrona) wire saw 40cm p/5	
Unit	box/5	Customs code (HS code)	90189084
GMDN	44987	GTIN	04046826088026
Risk category	1r (EU MDR)	UMDNS-Code	13-447
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;04;ai;59
Net weight (kg)	0,008	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media			
Complete picture		Detail images	
 			
size	40cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		