



## Product data sheet REF 57.70.02

11.12.2025

General data			
article number	57.70.02	Sørensen adult head rest only 340x200mm f adults	
Unit	piece	Customs code (HS code)	90189084
GMDN	35446	GTIN	04046826372064
Risk category	1 (EU MDD)	UMDNS-Code	11-958
Classification	1 (EU MDD)	Base UDI-DI	
Net weight (kg)	1,350	Operating instructions	G577000
Latex free	yes	Preparation instructions	MPK01
			
Medicon media	451.02.16 (brochure head rest system FiXiT)		
Complete picture		Detail images	
			
size	340x200mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.3535 (AIMG3)		