





## Product data sheet REF 57.70.16

11.12.2025

General data			
article number	57.70.16	two-spike support f adults w quick lock	
Unit	piece	Customs code (HS code)	90189084
GMDN	35446	GTIN	04046826372071
Risk category	1 (EU MDD)	UMDNS-Code	11-958
Classification	1 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,025	Operating instructions	G051570
Latex free	yes	Preparation instructions	MPK01
			
Medicon media	451.02.16 (brochure head rest system FiXiT)		
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.3535 (AIMG3)		