



Product data sheet REF 57.85.01

11.12.2025

| General data | | | |
|--|-----------------------------------|--------------------------|----------------|
| article number | 57.85.01 | Raney scalp clip p/10 | |
| Unit | box/10 | Customs code (HS code) | 90189084 |
| GMDN | 46953 | GTIN | 04046826088255 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 10-909 |
| Classification | 7 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,008 | Operating instructions | |
| Latex free | yes | Preparation instructions | MPK01 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4310 (X10CrNi18-8) | | |