



## Product data sheet REF 57.85.03A

11.12.2025

General data			
article number	57.85.03A	Raney scalp clip blue single-use p/20x10	
Unit	box/200	Customs code (HS code)	90189084
GMDN	46953	GTIN	04046826386894
Risk category	2a (EU MDD)	UMDNS-Code	10-909
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,220	Operating instructions	G578503A
Latex free	yes	Preparation instructions	
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 2px 5px;"><b>STERILE</b></div> <div style="border: 1px solid black; padding: 2px 5px;"><b>EO</b></div>   </div>			
Medicon media			
Complete picture		Detail images	
			
Storage conditions			
5 °C	 30 °C	10 %	 55 %
			
			
Materials			
	<b>Material number (DIN designation)</b>		
<b>Main part</b>	MAKROLON RX1805 (Polycarbonat)		