





Product data sheet REF 58.24.17

11.12.2025

General data			
article number	58.24.17	suct tube Ø 2.0mm 12cm conical mall	
Unit	piece	Customs code (HS code)	90189084
GMDN	38749	GTIN	04046826088606
Risk category	2a (EU MDD)	UMDNS-Code	10-566
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,025	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK03-1
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 2.0mm 12cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		