






## Product data sheet REF 58.24.20

11.12.2025

General data			
article number	58.24.20	suct tube Ø 4.0mm 12cm conical mall	
Unit	piece	Customs code (HS code)	90189084
GMDN	38749	GTIN	04046826088637
Risk category	2a (EU MDD)	UMDNS-Code	10-566
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,025	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK03-1
			
Medicon media			
Complete picture		Detail images	
 			
size	Ø 4.0mm 12cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4301 (X5CrNi18-10)		
<b>Handle</b>	1.4301 (X5CrNi18-10)		
<b>Others</b>	1.4305 (X8CrNiS18-9)		