







## Product data sheet REF 58.30.17

11.12.2025

General data			
article number	58.30.17	Yasargil raspatory 18.5cm 2.5mm slight cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	35559	GTIN	04046826088743
Risk category	1r (EU MDR)	UMDNS-Code	11-084
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;03;ai;4W
Net weight (kg)	0,016	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT ) 450.00.46 (catalogue sinuscopy )		
Complete picture		Detail images	
 			
size	18.5cm 2.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		