







Product data sheet REF 58.30.20

11.12.2025

General data			
article number	58.30.20	Jacobson buttoned probe ang 18.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	33433	GTIN	04046826088767
Risk category	1r (EU MDR)	UMDNS-Code	15-237
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;06;aj;6M
Net weight (kg)	0,033	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT) 451.02.51 (ENT basic sets Dr M W Winter Instr sets f ENT surgery)		
Complete picture		Detail images	
 			
size	18.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		