



## Product data sheet REF 58.30.65

11.12.2025

General data			
article number	58.30.65	Yasargil lig guide 18.5cm ang 90°	
Unit	piece	Customs code (HS code)	90189084
GMDN	32864	GTIN	04046826088927
Risk category	1r (EU MDR)	UMDNS-Code	12-336
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,052	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
			
size	18.5cm 90°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		