







Product data sheet REF 58.32.18

11.12.2025

General data			
article number	58.32.18	micro hook 1.5mm bl 90° 23cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	46768	GTIN	04046826089054
Risk category	1r (EU MDR)	UMDNS-Code	12-029
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,019	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
 			
size	1.5mm 90° 23cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		