







Product data sheet REF 58.32.64

11.12.2025

| General data | | | |
|--|-----------------------------------|--|-------------------------|
| article number | 58.32.64 | Rhoton micro hook 2mm 90° semi sh 19cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 46768 | GTIN | 04046826089238 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-029 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,020 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK03-1 |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 2mm 90° 19cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Handle | 1.4305 (X8CrNiS18-9) | | |