







## Product data sheet REF 58.47.65

11.12.2025

General data			
article number	58.47.65	Yasargil Stand Ti clip fenestr 5.0mm ang 5mm te	
Unit	piece	Customs code (HS code)	90219090
GMDN	34958	GTIN	04046826091439
Risk category	3 (EU MDD)	UMDNS-Code	10-905
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,003	Operating instructions	G585408
Latex free	yes	Preparation instructions	MPK18
  			
Medicon media			
Complete picture		Detail images	
			
size	5.0mm 5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		