







Product data sheet REF 58.55.55

11.12.2025

General data			
article number	58.55.55	Yasargil Stand clip 90° ang 5mm p	
Unit	piece	Customs code (HS code)	90219090
GMDN	34958	GTIN	04046826092597
Risk category	3 (EU MDD)	UMDNS-Code	10-905
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,005	Operating instructions	G585408
Latex free	yes	Preparation instructions	MPK18
  			
Medicon media			
Complete picture		Detail images	
			
size	90° 5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.4711 (CoCr20Ni15Mo)		