







## Product data sheet REF 58.56.92

11.12.2025

General data			
article number	58.56.92	Yasargil Stand Ti clip appl fcps rot 23cm swivelling	
Unit	piece	Customs code (HS code)	90189084
GMDN	32591	GTIN	04046826092979
Risk category	1r (EU MDR)	UMDNS-Code	10-904
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R005;01;ak;5L
Net weight (kg)	0,030	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media			
Complete picture		Detail images	
 			
size	23cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4021 (X20Cr13)		
<b>Screw</b>	1.4305 (X8CrNiS18-9)		
<b>Others</b>	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		