







Product data sheet REF 58.59.89

11.12.2025

| General data | | | |
|---|-------------------------------------|--|----------------|
| article number | 58.59.89 | Yasargil Stand Ti clip fenestr 5.0mm ang 7mm p | |
| Unit | piece | Customs code (HS code) | 90219090 |
| GMDN | 34958 | GTIN | 04046826094089 |
| Risk category | 3 (EU MDD) | UMDNS-Code | 10-905 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,003 | Operating instructions | G585408 |
| Latex free | yes | Preparation instructions | MPK18 |
|    | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 5.0mm 7mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) | | |