



Product data sheet REF 58.60.00

11.12.2025

General data			
article number	58.60.00	tray f Yasargil clips 66 cavities	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826094157
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	1,360	Operating instructions	G585408
Latex free		Preparation instructions	KEIN MP
			
Medicon media			
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	PPSU (Polyphenylsulfon) 1.4301 (X5CrNi18-10)		