







Product data sheet REF 58.61.20

11.12.2025

General data			
article number	58.61.20	Yasargil Stand Ti clip cvd 20mm p	
Unit	piece	Customs code (HS code)	90219090
GMDN	34958	GTIN	04046826094577
Risk category	3 (EU MDD)	UMDNS-Code	10-905
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,003	Operating instructions	G585408
Latex free	yes	Preparation instructions	MPK18
  			
Medicon media			
Complete picture		Detail images	
			
size	20mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		