







Product data sheet REF 58.65.07

11.12.2025

General data			
article number	58.65.07	Yasargil Stand Ti clip 90° ang 7mm te	
Unit	piece	Customs code (HS code)	90219090
GMDN	34958	GTIN	04046826095512
Risk category	3 (EU MDD)	UMDNS-Code	10-905
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,003	Operating instructions	G585408
Latex free	yes	Preparation instructions	MPK18
  			
Medicon media			
Complete picture		Detail images	
			
size	90° 7mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		