







Product data sheet REF 59.47.18

11.12.2025

General data			
article number	59.47.18	ti micro fcps Ø 1mm 18cm tips diamond dust coated	
Unit	piece	Customs code (HS code)	90189084
GMDN	62466	GTIN	04046826096274
Risk category	1r (EU MDR)	UMDNS-Code	15-621
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;03;ah;48
Net weight (kg)	0,028	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.53 (catalogue cardio) 450.00.60 (catalogue plastic surgery) 451.04.41 (brochure titanium instruments f Cardiac and Vascular Surgery)		
Complete picture		Detail images	
 			
size	Ø 1mm 18cm	Coating	Artificially produced diamond coating on working ends of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		