



Product data sheet REF 60.00.63

11.12.2025

General data			
article number	60.00.63	Bonn fixation ring 3/4 opened Ø 13mm	
Unit	piece	Customs code (HS code)	90185010
GMDN	60909	GTIN	04046826096618
Risk category	1 (EU MDR)	UMDNS-Code	17-933
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;04;XN
Net weight (kg)	0,040	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media	450.00.03 (general catalogue the art of surgery)
---------------	--

Complete picture

Detail images



size	Ø 13mm 11cm	Coating	
------	-------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)
Handle	1.4305 (X8CrNiS18-9)