







Product data sheet REF 60.22.02

11.12.2025

General data			
article number	60.22.02	Tyrell iris hook sh 12.5cm	
Unit	piece	Customs code (HS code)	90185010
GMDN	32754	GTIN	04046826097899
Risk category	1r (EU MDR)	UMDNS-Code	13-381
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,034	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
 			
size	12.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		