



Product data sheet REF 60.22.17

11.12.2025


General data			
article number	60.22.17	Sinsky lens positioning hook 11.5cm ang dw Ø 0.20mm	
Unit	piece	Customs code (HS code)	90185010
GMDN	32754	GTIN	04046826097967
Risk category	1r (EU MDR)	UMDNS-Code	13-381
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,025	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1



Medicon media	450.00.03 (general catalogue the art of surgery)
---------------	--

Complete picture

Detail images



size	11.5cm Ø 0.20mm	Coating	
------	-----------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)
Handle	1.4021 (X20Cr13)