








Product data sheet REF 60.32.02

11.12.2025

General data			
article number	60.32.02	Wecker iris spatula 2mm 13.5cm sterling silver	
Unit	piece	Customs code (HS code)	90185010
GMDN	39610	GTIN	04046826098353
Risk category	1r (EU MDR)	UMDNS-Code	13-648
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,012	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	2mm 13.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	SILBER_800 (AgCu_80-20) 1.4021 (X20Cr13)		