



Product data sheet REF 60.33.03

11.12.2025

General data			
article number	60.33.03	Anel tear duct cann 90° cvd 23 Gauge 4cm LL	
Unit	piece	Customs code (HS code)	90185010
GMDN	17899	GTIN	04046826098582
Risk category	1 (EU MDR)	UMDNS-Code	10-574
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;02;XG
Net weight (kg)	0,009	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 453.03.01 (Blepharoplasty - Instruments for blepharoplasty Brochure)		
Complete picture		Detail images	
 			
size	90° 23 Gauge 4cm	Coating	Electroplated chrome layer on the surface
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	2.0401 (CuZn39Pb3)		