








Product data sheet REF 60.33.23

11.12.2025

General data			
article number	60.33.23	Knolle pearce lens loop w irrig 23 Gauge 5cm	
Unit	piece	Customs code (HS code)	90185010
GMDN	12319	GTIN	04046826098681
Risk category	1r (EU MDR)	UMDNS-Code	13-381
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;02;ai;4K
Net weight (kg)	0,015	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	23 Gauge 5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		