



Product data sheet REF 60.73.50

11.12.2025

| General data | | | |
|--|--|---------------------------|-------------------------|
| article number | 60.73.50 | Barraquer iris scrs 5.5cm | |
| Unit | piece | Customs code (HS code) | 90185010 |
| GMDN | 13488 | GTIN | 04046826100896 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-488 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;01;ag;2U |
| Net weight (kg) | 0,026 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 5.5cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Screw | 1.4034 (X46Cr13) 1.4305 (X8CrNiS18-9) | | |