







## Product data sheet REF 62.16.03

11.12.2025

General data			
article number	62.16.03	Wagener ear hook 14cm 3.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	13039	GTIN	04046826102944
Risk category	1r (EU MDR)	UMDNS-Code	15-258
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,016	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
 			
size	14cm 3.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		