





Product data sheet REF 62.21.48

11.12.2025

| General data | | | |
|---|-----------------------------------|--------------------------|----------------|
| article number | 62.21.48 | suct adapter luer tip | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 38749 | GTIN | 04046826103729 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 10-571 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,013 | Operating instructions | |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | 450.00.37 (catalogue ENT) | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4305 (X8CrNiS18-9) | | |