





## Product data sheet REF 62.22.35

11.12.2025

| General data   |                                   |                           |                |
|--|-----------------------------------|---------------------------|----------------|
| article number   | 62.22.35                          | Zöllner suct cann Ø 0.5mm |                |
| Unit   | piece                             | Customs code (HS code)    | 90189084       |
| GMDN   | 38749                             | GTIN                      | 04046826103781 |
| Risk category  | 2a (EU MDD)                       | UMDNS-Code                | 10-571         |
| Classification   | 6 (EU MDD)                        | Base UDI-DI               |                |
| Net weight (kg)  | 0,001                             | Operating instructions    |                |
| Latex free   | yes                               | Preparation instructions  | MPK02          |
|   |                                   |                           |                |
| Medicon media  | 450.00.37 (catalogue ENT )        |                           |                |
| Complete picture   |                                   | Detail images             |                |
|  |                                   |                           |                |
| size   | Ø 0.5mm                           | Coating                   |                |
| Storage conditions   |                                   |                           |                |
| No special requirements  |                                   |                           |                |
| Materials  |                                   |                           |                |
|  | Material number (DIN designation) |                           |                |
| Main part  | 1.4305 (X8CrNiS18-9)              |                           |                |
| Others   | 1.4301 (X5CrNi18-10)              |                           |                |