





## Product data sheet REF 62.23.30

11.12.2025

General data			
article number	62.23.30	suct and irrig tube Ø 3.0+2.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	38749	GTIN	04046826103897
Risk category	2a (EU MDD)	UMDNS-Code	10-571
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,017	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.37 (catalogue ENT ) 450.00.46 (catalogue sinuscscopy )		
Complete picture		Detail images	
			
size	Ø 3.0+2.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		