








Product data sheet REF 62.46.11

11.12.2025

General data			
article number	62.46.11	Schuknecht wire crimper 7.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32886	GTIN	04046826105006
Risk category	1r (EU MDR)	UMDNS-Code	15-221
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R005;03;ak;6A
Net weight (kg)	0,033	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
  			
size	7.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		