



Product data sheet REF 62.60.70

11.12.2025

General data			
article number	62.60.70	Lempert elevator heavy 18cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35559	GTIN	04046826105310
Risk category	1r (EU MDR)	UMDNS-Code	15-211
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;03;ai;4W
Net weight (kg)	0,038	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1



Medicon media 450.00.37 (catalogue ENT)

Complete picture

Detail images



size 18cm Coating

Storage conditions

No special requirements

Materials

Material number (DIN designation)

Main part
Handle 1.4021 (X20Cr13)
1.4301 (X5CrNi18-10)