







## Product data sheet REF 63.06.70

11.12.2025

General data			
article number	63.06.70	hook 15cm 90° 2.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	13039	GTIN	04046826105808
Risk category	1r (EU MDR)	UMDNS-Code	15-258
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,034	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
 			
size	15cm 90° 2.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		