







Product data sheet REF 63.10.03

11.12.2025

| General data | | | |
|--|--|-------------------------------------|-------------------------|
| article number | 63.10.03 | Guilford Wright flap knife 15cm 2mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 12245 | GTIN | 04046826106461 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-245 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;02;ag;37 |
| Net weight (kg) | 0,034 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK01-1 |
|   | | | |
| Medicon media | 450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 15cm 2mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |