








Product data sheet REF 63.11.03

11.12.2025

General data			
article number	63.11.03	Goldmann elevator 15cm 2.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35559	GTIN	04046826106591
Risk category	1r (EU MDR)	UMDNS-Code	12-245
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,034	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
  			
size	15cm 2.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		