







Product data sheet REF 63.14.10

11.12.2025

General data			
article number	63.14.10	Hough depth gauge 17cm grad 3.5/4.0/4.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32845	GTIN	04046826106997
Risk category	1r (EU MDR)	UMDNS-Code	15-237
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;06;aj;6M
Net weight (kg)	0,033	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
 			
size	17cm 3.5/4.0/4.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		