






Product data sheet REF 63.90.27

11.12.2025

General data			
article number	63.90.27	trephine stainless Ø 2.7mm 44mm	
Unit	piece	Customs code (HS code)	90184910
GMDN	14147	GTIN	04046826109028
Risk category	2a (EU MDD)	UMDNS-Code	10-520
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,002	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
 			
size	Ø 2.7mm 44mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4197 (X20CrNiMoS13-1)		