





Product data sheet REF 64.54.15

11.12.2025

| General data | | | |
|---|-----------------------------------|--------------------------|--|
| article number | 64.54.15 | diamond burr 95mm 1.4mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 46392 | GTIN | 04046826109448 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 15-883 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,004 | Operating instructions | |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | 450.00.46 (catalogue sinuscopy) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 95mm 1.4mm | Coating | Artificially produced diamond coating on working ends of instruments |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4112 (X90CrMoV18) | | |